

Research Article

Colistin broth disk elution for the detection of Colistin resistance among Carbapenem-resistant *Klebsiella pneumoniae*

Elución en caldo con disco de colistina para la detección de resistencia a la colistina entre *Klebsiella pneumoniae* resistente a carbapenémicos

Phu Truong Thien^{1,2} https://orcid.org/0000-0003-0868-2811 Ngoc Tran Bich¹* https://orcid.org/0000-0002-2255-6398 Hoang Pham Huy¹ https://orcid.org/0009-0008-2090-5268 Van Tran Thi Hue¹ https://orcid.org/0009-0007-0677-2553 Hoang Phan Thanh¹ https://orcid.org/0009-0003-6737-7216 Binh Pham Thai¹ https://orcid.org/0009-0002-2475-8006 Lam Tran Thi Ngoc¹ https://orcid.org/0009-0005-9125-4771 Truc Nguyen Thi Thanh¹ https://orcid.org/0009-0008-6117-4270 Mai Le Phuong² https://orcid.org/0009-0001-9625-2522 Dat Ngo Quoc¹ https://orcid.org/0000-0003-1461-0216

¹University of Medicine and Pharmacy at Ho Chi Minh City. Ho Chi Minh City, Vietnam. ²Cho Ray Hospital. Ho Chi Minh City, Vietnam.

*Author for correspondence. Email: bichngoctran@ump.edu.vn

ABSTRACT

Introduction: The emergence of multidrug-resistant pathogens, especially carbapenem-resistant *Klebsiella pneumoniae* (CRKP), is a pressing public health issue. The limited therapeutic choices further compound this. The escalating resistance to colistin, a crucial last-resort antibiotic, underscores the need for precise and dependable antimicrobial susceptibility testing (AST) methods.

Objectives: To evaluate the performance of the Colistin Broth Disk Elution (CBDE) method in comparison to the reference Broth Microdilution (rBMD) method for detecting colistin resistance in CRKP isolates.

Methods: A cross-sectional study involving 129 non-duplicate CRKP isolates was conducted at Cho Ray Hospital from October 2023 to December 2023. Colistin susceptibility was assessed using both the BMD and CBDE methods. The agreement between the two methods was measured in terms of categorical agreement (CA) and essential agreement (EA) based on Clinical and Laboratory Standards Institute (CLSI) criteria.

Results: Of the 129 CRKP isolates, 23.26% were identified as colistin-resistant by both BMD and CBDE methods. The CA between the two methods was 100%, and the EA was 91.47%, with no major errors (ME) or very major errors (VME) detected.

Conclusion: The CBDE method has shown a high level of agreement with the reference BMD method, making it a reliable alternative for the routine detection of colistin resistance in microbiology laboratories. The study's findings strongly support using CBDE as a practical and effective AST method for colistin, thereby contributing to antibiotic stewardship.

Keywords: antimicrobial susceptibility testing; carbapenem resistance; colistin; *Klebsiella pneumoniae*; Microdilution.

RESUMEN

Introducción: El aumento de patógenos multirresistentes, especialmente *Klebsiella pneumoniae* resistente a carbapenémicos (CRKP), representa una preocupación para la salud pública, agravada por opciones de tratamiento limitadas. La creciente resistencia a la colistina, un antibiótico de último recurso, resalta la necesidad de métodos precisos y confiables de pruebas de susceptibilidad antimicrobiana (AST).

Objetivos: Evaluar la eficacia del método de elución en caldo con disco de colistina (CBDE) en comparación con el método estándar de microdilución en caldo (BMD) para detectar resistencia a la colistina en aislados de CRKP.

Métodos: Estudio transversal con 129 aislados únicos de CRKP, en el Hospital Cho Ray, desde octubre a diciembre de 2023. La susceptibilidad a la colistina se evaluó utilizando, tanto los métodos BMD como CBDE. La concordancia entre métodos se midió en términos de acuerdo categórico (CA) y acuerdo esencial (EA), basados en los estándares del Instituto de Estándares Clínicos y de Laboratorio (CLSI).

Resultados: Entre los 129 aislados de CRKP, el 23,26% se identificaron como resistentes a la colistina por ambos métodos. La CA entre BMD y CBDE fue del 100%, con un EA del 91,47%, y no se observaron errores mayores (ME) ni muy graves (VME).

Conclusión: El método CBDE demostró una alta concordancia con el BMD de referencia, se establece como una alternativa confiable para la detección rutinaria de resistencia a colistina en laboratorios de microbiología. Estos hallazgos respaldan el uso de CBDE como un método AST práctico y efectivo para la colistina y apoya los esfuerzos de gestión de antibióticos.

Palabras clave: colistina; *Klebsiella pneumoniae*; microdilución en caldo; prueba de susceptibilidad antimicrobiana; resistencia a carbapenémicos.

Received: 30/08/2024 Approved: 25/11/2024

INTRODUCTION

Multidrug resistance is a crucial problem threatening public health worldwide, increasing healthcare costs in diagnostic and therapy. Recently, Enterobacterales have received particular attention because of their high resistance to carbapenem and polymyxin B.⁽¹⁾ Among these bacteria, *K. pneumoniae* is a Gramnegative bacterium commonly found in healthcare settings and is responsible for many infections, such as pneumonia, bloodstream infections, and urinary tract infections.⁽²⁾ The ability of *K. pneumoniae* to rapidly acquire and disseminate resistance genes has created multidrug-resistant strains, causing difficulty in treating diseases. Carbapenem-resistant *K. pneumoniae* (CRKP) has been particularly problematic, often associated with high mortality rates and limited treatment options.⁽²⁾

http://scielo.sld.cu https://revmedmilitar.sld.cu

Bajo licencia Creative Commons (6)

The spreading of CRKP and carbapenem-resistant Enterobacterales led to the revival of colistin as a lastline treatment, even though its use requires extreme caution.⁽³⁾ Unfortunately, the use of colistin in both human and veterinary medicine has caused colistin resistance.^(4,5) Hence, it becomes necessary that bacteria be identified with colistin resistance in laboratories before considering using it for treatment.^(3,4) There is an increasing need for a standardized antimicrobial susceptibility test (AST) for colistin, vital for patient care and surveillance.

The only reference method for colistin AST recognized by the Clinical and Laboratory Standards Institute (CLSI) and European Committee on Antimicrobial Susceptibility Testing (EUCAST) is Broth Microdilution (BMD).^(6,7) Although BMD is a standard method, it is challenging to do routinely because of its laborious nature. However, CBDE, which has also assessed alternative approaches by CLSI, seems more suitable for routine application in microbiology laboratories due to its relatively simple materials and media.^(7,8) In light of this, this study aims to evaluate the performance of the CBDE method compared to the rBMD method for detecting colistin resistance in CRKP isolates.


METHOD

Study design and participants


A cross-sectional study was conducted at Cho Ray Hospital from October 2023 to December 2023, with 130 strains of *K. pneumoniae* collected from routine clinical specimens of blood, urine, sputum, pus, and body fluid. The figure 1 display the flow chart of participants throughout the study duration. These strains resisted at least one carbapenem antibiotic and had colistin-BMD results. Bacterial strains isolated for the second time in the same patient or determined to be from an outbreak, as notified by the infection control department, were excluded from the study. One bacterial strain was excluded due to contamination; there were 129 eligible participant strains.

Revista Cubana de Medicina Militar

rBMD



CBDE



Index test n= 129

Fig 1 - Flow chart of participants.

Data collection

Bacterial strains were identified using the Vitek MS System (BioMérieux, France), and carbapenemresistant strains were detected with the Vitek-2 Compact System (BioMérieux, France). After assessment, each isolate was preserved in Tryptic Soy Broth (TSB) with 15% glycerol for future use. Results were recorded and analyzed. The reference BMD was conducted according to the SensitiveTM FRCOL protocol. Cation-adjusted Mueller-Hinton broth procured from OxoidTM was used for the CBDE method, performed according to the CLSI M100 S34. Colistin sulfate salt was obtained from OxoidTM. *Pseudomonas aeruginosa* ATCC[®] 27853 and *Escherichia coli* ATCC[®] 25922 were used for quality control.

According to CLSI M100 S34, colistin is considered resistant to bacteria strains when the minimum inhibitory concentration (MIC) is $\geq 4 \ \mu g/mL$ and intermediate to bacteria strains when MIC is $\leq 2 \ \mu g/mL$.

Statistical Analysis

The MIC results of CBDE were compared with the BMD method. The reference BMD test and CBDE agreement were measured in terms of categorical agreement (CA) and essential agreement (EA). CA was

determined as the percentage of samples with results in the same category as the reference method, using all tested isolates as the denominator (n= 129). EA was calculated as the percentage of isolates with MIC values within \pm 1 log2 dilution or \pm 1 twofold dilution of the reference standard. According to CLSI M52 S15, any test with both CA and EA exceeding 90% is considered a reliable alternative to the reference test.⁽⁹⁾ Errors were classified as very major errors (VME: false-intermediate) and major errors (ME: false-resistant) by CBDE, with VMEs and MEs of \geq 3% deemed unacceptable.

Ethical considerations

The study collected frozen bacterial strains that met the selection criteria without patient intervention. The Ethics Committee of the University of Medicine and Pharmacy in Ho Chi Minh City approved this research under contract number 123/HĐĐĐ-ĐHYD.

RESULTS

Characteristics of the study participants

Table 1 shows that among the 129 CRKP isolates tested in the study, the respiratory tract remains the primary route of infection; 35.7% of CRKP was isolated from phlegm specimens. Most bacterial strains were identified as resistant to ertapenem, imipenem, and meropenem, accounting for 98.4%.



Characteristics	n	%		
Specimens				
Sputum	46	35.7		
Body fluid	34	26.4		
Blood	24	18.6		
Urine	15	11.6		
Pus	10	7.8		
Resistance to carbapenems				
Ertapenem	1	0.8		
Imipenem & Ertapenem	1	0.8		
Imipenem, Ertapenem & Meropenem	127	98.4		
Total	129	100		

 Table 1 - Characteristics of the study population

Prevalence of colistin resistance

Table 2 shows the distribution of MIC values in participant strains. The CBDE method, as the BMD method does, cannot determine exactly colistin MIC values at $< 1 \ \mu g/mL$ and $> 4 \ \mu g/mL$ levels. The first data row represents the MIC breakpoint as per CLSI standards. Of CBDE results, 11 isolates had MIC values exceeding ± 1 two-fold dilution of the reference standard (Fig. 2).

Table 2 - MIC distribution obtained by the reference BMD method and CBDE method

Minimum inhibitory concentration	CBDE MIC (µg/mL)			
rBMD MIC (µg/mL)	< 1	1	2	≥4
< 1	-	56	11	-
1	-	26	3	-
2	-	1	2	-
\geq 4	-	-	-	30

rBMD: The reference Broth Microdilution; CDBE: Colistin Disk Broth Elution method.





Fig 2 - A result of CBDE.

Table 3 shows that this study analyzed 129 CRKP isolates. Of these, 30 out of 129 participant strains were identified as colistin-resistant by the BMD and the CBDE methods, accounting for 23.26%. No false negative or false positive was detected.

Table 3 - Results of BMD and	CBDE for 129 CRKP isolates
------------------------------	----------------------------

Resistance	B	MD	CBDE	
Resistance	n	%	n	%
Intermediate (MIC $\leq 2 \ \mu g/mL$)	99	76.74	99	76.74
Resistant (MIC \ge 4 µg/mL)	30	23.26	30	23.26
Total	129	100	129	100

BMD: Broth Microdilution; CDBE: Colistin Disk Broth Elution method; MIC: Minimum inhibitory concentration.

The agreement of CBDE to reference BMD for detecting colistin resistance

Table 4 compares the CBDE results with the MIC values obtained by the BMD method. The CBDE method's CA and EA were within acceptable limits according to CLSI evaluation standards: CA was

100%, and EA was 91.47%. No ME or VME was detected in this study. The results suggested that the CBDE method was highly agreeable with the BMD method in detecting colistin resistance.

Method	CRKP (n = 129)		No. (%) of isolates achieving				CLSI acceptable limits	
	Ι	R	СА	EA	ME	VME	CA, EA	ME, VME
BMD	99	30	129/129	118/129	0 (0) 0 (0)	0 (0)	> 90%	< 3%
CBDE	99	30	(100)	(91,47)				

 Table 4 - Performance of CBDE in comparison to reference BMD method

CA: categorical agreement; EA: essential agreement; VME: very major errors (false-intermediate); ME: major errors (false-resistant).

DISCUSSION

The rise of carbapenem and colistin-resistant pathogens in clinical settings, coupled with the challenges associated with colistin susceptibility testing, is leading to a critical situation in treating severe infections. Implementing the joint CLSI-EUCAST recommended rBMD is challenging in diagnostic clinical microbiology laboratories, highlighting the need to evaluate alternative methods for therapeutic decision-making.⁽⁸⁾

A total of 129 CRKP, isolated from various clinical specimens, were included in this study. Colistin susceptibility was determined using BMD and CBDE methods based on CLSI M52 S15 recommendations.⁽⁹⁾ The CBDE method was compared with the reference BMD method based on evaluation criteria recommended by CLSI to assess the concordance between the two methods. The evaluation criteria included CA% and EA% \geq 90%, and ME% and VME% < 3%.

In this study, 23.26% (30/129) of the CRKP isolates were resistant to colistin, as determined by the reference BMD method. The CBDE method also identified 23.26% (30/129) of the isolates as colistin-resistant. Based on the collected data, the CBDE method demonstrated a CA% of 100% and an EA% of 91.47%, with no ME or VME recorded. These results indicate that the CBDE method possesses high sensitivity and specificity for determining colistin susceptibility in CRKP and CRE. It is an effective and straightforward method that could be considered a reliable alternative to the BMD test.

Recent studies worldwide have shown results that align well with our findings on the comparability of the BMD and CBDE methods. For instance, a multicenter study was conducted in 2019 to evaluate the efficacy of various testing methods for detecting colistin resistance, as outlined in the CLSI report by *Humphries* et al.⁽¹⁰⁾ In this study, the CBDE method achieved a CA of 96.8% for 348 Enterobacterales isolates, with VME was 2.5%, while no ME was detected. In comparison, a study by *Sujatha* et al.⁽¹¹⁾ in India reported a CA of 99% for the CBDE method, with an ME rate of only 1% compared to the reference BMD method. Similarly, *Kansak N* et al.⁽⁸⁾ in Turkey observed that the CBDE method achieved a perfect CA of 100%, and, consistent with our findings, no ME or VME was observed. The researchers concluded that the CBDE method demonstrates robust performance compared to the rBMD method, deeming it suitable for routine use in detecting colistin resistance. Furthermore, in a study conducted in India in 2023, *Rout* et al.⁽¹²⁾ reported values of 97.6% for CA, 94.7% for EA, 2.26% for ME, and 1.9% for VME in the CBDE method.

In conclusion, the CBDE method, as implemented in the current study and corroborated by other studies globally, meets the CLSI criteria for comparability with the reference BMD method. The current study provides additional scientific evidence supporting the efficacy of the CBDE method, positioning it as a viable alternative to the BMD test in laboratories of various scales in Vietnam.

There are also some limitations: this research was about the colistin phenotype and showed no mobile cases of colistin resistance. Furthermore, the sample size should come from a multicenter to be more significant in detecting colistin resistance, especially in *E. coli*.

The CBDE method has shown a high level of agreement with the reference BMD method, making it a reliable alternative for the routine detection of colistin resistance in microbiology laboratories. The study's findings strongly support using CBDE as a practical and effective AST method for colistin, thereby contributing to antibiotic stewardship.

Acknowledgments

We especially thank University of Medicine and Pharmacy at Ho Chi Minh City for funding and Cho Ray Hospital for providing the facilities for this study. Furthermore, we are honored to aid samples from the Microbiology Department, Cho Ray Hospital, in implementing our study.

BIBLIOGRAPHIC REFERENCES

1. Muray CJ, Ikuta KS, Sharara S, Swetschinski L, Aguilar GR, Gray A, et al. Global burden of bacterial antimicrobial resistance in 2019: a systematic analysis [Internet]. The Lancet. 2022; 399(10325):629-55. DOI: 10.1016/S0140-6736(21)02724-0

Grundmann H, Glasner C, Albiger B, Aanensen DM, Tomlinson CT, Andrasevic AT, et al.
 Occurrence of carbapenemase-producing Klebsiella pneumoniae and Escherichia coli in the European survey of carbapenemase-producing Enterobacteriaceae (EuSCAPE): a prospective, multinational study [Internet]. The Lancet Infectious Diseases. 2017; 17(2):153–63. DOI: 10.1016/S1473-3099(16)30257-2

3. Bir R, Gautam H, Arif N, Chakravarti P, Verma J, Banerjee S, et al. Analysis of colistin resistance in carbapenem-resistant Enterobacterales and XDR Klebsiella pneumoniae [Internet]. Ther Adv Infect Dis. 2022; 26(9):1-15. DOI: 10.1177/20499361221080650

4. Zhou K, Luo Q, Wang Q, Huang C, Lu H, Rossen JW, el al. Silent transmission of an IS1294bdeactivated mcr-1 gene with inducible colistin resistance [Internet]. International Journal of Antimicrobial Agents. 2018; 51(6):822-8. DOI: 10.1016/j.ijantimicag.2018.01.004

 Poirel L, Jayol A, Nordmann P. Polymyxins: Antibacterial Activity, Susceptibility Testing, and Resistance Mechanisms Encoded by Plasmids or Chromosomes [Internet]. Clin Microbiol Rev. 2017; 30(2):557-96. DOI: 10.1128/cmr.00064-16

6. Clinical and Laboratory Standards Institute (CLSI). Performance Standards for Antimicrobial Susceptibility Testing 34th ed. CLSI supplement M100. Wayne, PA: CLSI; 2024.

7. European Committee on Antimicrobial Susceptibility Testing (EUCAST). Breakpoint tables for interpretation of MICs and zone diameters [Internet]. Version 11.0. Basel: EUCAST; 2021 [access:



28/02/2024]. Available from:

https://www.eucast.org/fileadmin/src/media/PDFs/EUCAST_files/Breakpoint_tables/v_11.0_Breakpoint_tables.pdf

 Kansak N, Arici N, Uzunoner Y, Adaleti R, Aksaray S, Gonullu N. Evaluation of Broth Disk Elution Method to Determine Colistin Resistance in Klebsiella pneumoniae and Escherichia coli Strains [Internet]. Clinical Laboratory. 2023; 69(2): 221008. DOI: 10.7754/clin.lab.2022.221008
 Clinical and Laboratory Standards Institute (CLSI). Verification of Commercial Microbial Identification and Antimicrobial Susceptibility Testing Systems [Internet]. 15th ed. CLSI supplement M52. Wayne, PA: CLSI; 2015.

 Humphries RM, Green DA, Schuetz AN, Bergman Y, Lewis S, Yee R, et al. Multicenter evaluation of Colistin Broth Disk Elution and colistin agar test: a report from the Clinical and Laboratory Standards Institute [Internet]. Journal of Clinical Microbiology. 2019; 57(11):10-1128. DOI: 10.1128/jcm.01269-19

11. Sujatha SR, Rajasekhar D, Badveti S, Manthravadi KK. Evaluation of Vitek2 and Colistin broth disk elution test in comparison with Micro broth dilution for Susceptibility testing of Colistin among Carbapenem Resistant Enterobacterales [Internet]. Indian J Public Health Res Dev. 2022; 13(4):290-5. DOI: 10.37506/ijphrd.v14i4.18632

12. Rout B, Dash SK, Kumar Sahu K, Behera B, Praharaj I, Otta S. Evaluation of different methods for in vitro susceptibility testing of colistin in carbapenem-resistant Gram-negative bacilli [Internet]. Access Microbiol. 2023; 5(10): 000595. DOI: 10.1099/acmi.0.000595.v3

Conflict of interest

The authors declare that they have no competing interests.

Financial information

This research was funded by University of Medicine and Pharmacy at Ho Chi Minh City under contract number 69/2024 HĐ-ĐHYD dated 15/04/2024.



Authors' Contributions

Conceptualization: Ngoc Tran Bich, Phu Truong Thien. Data curation: Hoang Pham Huy, Ngoc Tran Bich, Phu Truong Thien. Formal analysis: Ngoc Tran Bich, Hoang Pham Huy, Lam Tran Thi Ngoc, Truc Nguyen Thi Thanh, Phu Truong Thien. Research: Hoang Pham Huy, Hoang Phan Thanh, Ngoc Tran Bich, Phu Truong Thien. Methodology: Hoang Pham Huy, Hoang Phan Thanh, Ngoc Tran Bich, Phu Truong Thien. Supervision: Phu Truong Thien, Dat Ngo Quoc, Van Tran Thi Hue. Validation: Phu Truong Thien, Ngoc Tran Bich , Binh Pham Thai. Writing – original draft: Ngoc Tran Bich, Phu Truong Thien.

Data availability

This research data is confidential according to the applicable confidentiality agreements and regulations and, therefore, cannot be publicly displayed or shared. The data are securely stored at the Integrated Planning Department at Cho Ray Hospital. Access to these data requires proper authorization. If you have any questions or need further information, please contact Ngoc Tran Bich at <u>bichngoctran@ump.edu.vn</u>.